

Smile Makers Dental Laboratory

(Lab) Removable

(561) 966-8262

2540 Metrocentre Blvd, STE #4, West
Palm Beach, FL 33407

smilemakers@hizir3d.com

www.smilemakersfl.com

Internal Use Only

Rx

Pan # _____

Due Date: _____

Lab

Lab: _____
Address: _____
Phone: _____ Date: _____

PATIENT

First Name*: _____
Last Name*: _____
Appointment: _____

CASE ENCLOSURE

☐ Digital Scan ☐ Model ☐ Pics ☐ Impression ☐ Parts

Other: _____

CASE TYPE

☐ Metal Frame (CoCr) ☐ Metal Frame (Ti) ☐ Peek ☐ Vitallium

Kennedy Class: I II III IV Mod: _____ Edentulous Span: _____

Major Connector*:

Palatal Plate (Full/Partial) | A-P Bar | Single Palatal | Horseshoe | Lingual Bar | Lingual Plate | Double Lingual | Labial Bar

Clasp Types*:

1.Akers 2.RPI 3.RPA 4.Ring 5.Ball 6.Wrought Wire 7.Back Action 8.Reverse 9.Half&Half 10.Embrasure 11.Combination
12.Butterfly 13.Others: _____

DENTAL CHART

· Universal Numbering · Mark each tooth: **M** Missing **R** Rest Seat

Upper #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
M / R																
Lower #	32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17
M / R																

NOTE:

Signature: _____

Dental Lab License Number: _____

OUR STANDARD

Intaglio Surface smooth, well-polished standard — free of sharp edges and surface roughness for maximum patient comfort.

External Stippling light, low-profile finish. Deep or prominent texture is not used, minimizing soft-tissue irritation and patient complaints.

Framework Dimensions optimal balance of structural strength and patient comfort — never undersized, never bulky.

Clasp Arm Design adequate length and flexibility to deliver reliable retention without concentrating stress on the abutment tooth.

Opposing Components opposing components are positioned to clear the opposing dentition through all excursive movements.

Clasp Height & Position placed within the undercut — below the survey line — for dependable retention and proper soft-tissue clearance.